# MTC

Original

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D  
E

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Instructions

E

We are pleased that you decided to buy a highly developed piece of equipment from SCHICK and would like to wish you every success when working with your new MTC-control unit.

We wrote these instructions to enable you to get accustomed to your new piece of equipment and to provide you with the correct operating and maintenance instructions.

| Index  | Page |
|--|------|
| 1. Safety Instructions.....                      | 3    |
| Caution!   |      |
| 2. Range of applications.....                    | 3    |
| 3. Summary/components control unit cpl.....      | 4    |
| 4. Accessories/motorhandpieces and spindles..... | 4    |
| 5. Initial start up/Installation/Operation.....  | 5    |
| 6. Errors.....                                   | 6    |
| 7. Technical Data.....                           | 6    |
| 8. Declaration of Conformity.....                | 7    |

# 1. Safety Instructions

- 1.1 Ascertain that your mains supply coincides with the data on the rating plate.
- 1.2 MTC units are not suitable for the following applications:
  - in areas where there is a risk of explosion
  - for medical applications
  - for working on moist materials
- 1.3 Ensure that all regulatory requirements are observed during use
  - always wear protective glasses
  - provide enough light at the working place
  - use dust suction
- 1.4 Under no circumstances should the motor handpiece be cleaned with compressed air.

## Caution!

- When using rotary instruments, do not exceed the maximum speeds laid down by their manufacturer.
- Repairs and other technical procedures must only be carried out by suitable qualified personnel, authorized by SCHICK.
- SCHICK do not guarantee the MTC unit should it not have been used in accordance with the instructions.

# 2. Range of applications

The MTC control units in connection with SCHICK motors are designed for universal use in Industrial and handicraft business. It is designed for manual and hand-held working of various materials for example deburring, grinding, drilling and milling.

With its speed range from 1.000 - 50.000 rpm and high torque the SCHICK MTC control unit allows to work where strength, endurance and precision is required.

For the MTC units there are many motorhandpieces and spindles.

The powerful MTC microprocessor allows - through optimal motor control - fatigue-free work and making it economical to use over long periods. The progressive processor controller constantly monitors the unit for overloading. A second thermal fuse, which is independent of the computer, provides further protection. Thus the MTC unit offers a maximum of safety.

Conditions of environment:

- interior 5° - 40° C.
- up to 2,000 meter over sea level

Categorie of overvoltage: II

Grade of pollution: 2

WEEE-Rg.-Nr. DE 78620387

### 3. Summary Components control unit cpl. art.-no. 10280

MTC control unit art.no.: 10240


mains cable art.no.: 2160

foot switch art. no.: 6370 optional

### 4. Accessories: spindles and motorhandpieces

W 90 art.no.: 10101  
contra angle handpiece 90°


W 50 art.no.: 10103  
contra angle handpiece 50°

MT W Profi art.no.: 10105  
motor  
(with handpiece rack and tools)

MT Profi art.no.: 10210  
motorhandpiece MT Profi for unit MTC  
and MTQ (with handpiece rack and tools)

MT 6 - 1 art.no.: 10295  
spindle 6 mm for MTC and MTQ  
(incl. tools)

MT 6 art.no.: 10170  
motorhandpiece 6 mm for unit MTC  
and MTQ (with handpiece rack and tools)

IS 6 art.no.: 10200  
motor spindle 6 mm for  
unit MTC and MTQ  
(incl. tools)

IS 6 - 1 art.no.: 10200/1  
motor spindle 6 mm  
suitable for cutting fluid for  
MTC and MTQ  
(incl. tools)


BS I art.no.: 10290  
drilling spindle  
(incl. tools)

## 5. Initial start up / Installation

### 5.1 Initial start up

Connect the motorhandpiece to the motorhandpiece socket (1) (bayonet lock).  
Connect the foot switch to the foot switch socket (2). Plug the mains cable into the power supply socket (3) on the control unit.

 Basically, first attach the handpiece with the control unit before connecting the mains cable and switching the control mechanism on.

### 5.2 Operation:

The motorhandpiece can be started with the optional foot switch or directly through the speed controller.

1. Turn the mains switch to „on“
2. Start the motorhandpiece with the foot switch  
Select desired speed using the speed control on the control unit. Press and hold the foot switch.  
The motor stops when the foot switch is released
3. Start the motorhandpiece without the foot switch  
Set speed control to zero. The handpiece starts when the speed controller is started from zero.

The digital speed display shows the preset speed when the motorhandpiece stopped (see the red dot in front of the first digit) when the motor runs it shows the current speed (there is no red dot).

## 6. Errors

If the handpiece is overloaded or blocked, the control unit switches off automatically to prevent further damages.

Make sure that the tool can rotate freely in the handpiece. After releasing the foot controller or speed controller to neutral position the unit can continue operation.

If an error cannot be removed following the above description, please contact an authorised service partner or Schick directly.

## 7. Technical Data

dimensions MTC control unit:

|         |         |
|---------|---------|
| width:  | 225 mm  |
| height: | 110 mm  |
| depth:  | 170 mm  |
| weight: | 2,9 kgs |

|                    |              |
|--------------------|--------------|
| operating voltage: | AC 230 Volts |
| Frequency:         | 50/60 Hz     |
| output:            | 230 Watt     |

|                   |  |
|-------------------|--|
| speed range :     | 1.000 - 50.000 min                                 |
| torque:           | max: 9 Ncm (when using the optional Schick motors) |
| protection class: | I  |

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declare herewith, that the products

MTC - control unit cpl. art.-no. 10280 consisting of  
MTC - control unit art.-no. 10240

optionally with

foot switch on/off art.-no. 6370  
contra angle handpiece 90° art.-no. 10101  
contra angle handpiece 50° art.-no. 10103  
MT W Profi motor art.-no. 10105  
MT Profi motorhandpiece art.-no. 10210  
MT 6-1 industry spindle art.-no. 10295  
MT 6 motorhandpiece art.-no. 10170  
IS 6 industry spindle art.-no. 10200  
IS 6-1 industry spindle art.-no. 10200/1  
BS I drilling spindel art.-no. 10290

are in conformity with the following provisions of Directive:

2006/42/EG (machinery directive)  
2014/30/EU (EMC directive)  
2011/65/EU (RoHS)

Name and address of  
person in charge:

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Schemmerhofen, March 2017

W. Schick  
manager

Subject to technical modifications

This unit complies with the current VDE (German association of electrical technicians)  
regulations concerning safety and suppression.

These instructions should be readily accessible and are best kept close  
to the unit itself.

We would like to take this opportunity to advise you that a proper repair service and  
suitable qualified personnel are required for such highly developed technical equipment.  
SCHICK guarantees to carry out perfect repairs using original spare parts.

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Schick GmbH